

SOURCE CODE

1. Hardware Firmware (Arduino/C++)

File Name: MedLink_Firmware.ino

Description: This code runs on the Arduino Uno/ESP32. It controls the servo motor for pill dispensing, reads the Pulse Sensor for vitals, and communicates with the MedLink IoT app.

C++

```
/*
 * Project Topic: IoT Based Smart Pillbox & Vital Health Monitoring
 * App Name: MedLink IoT
 * Created By: Group 51-55 (TCET)
 * * Description:
 * This firmware controls the servo motor for pill dispensing, reads
 * real-time health data from the MAX30102 sensor, and updates the
 * OLED display. It communicates with the cloud via ESP32.
 */

#include <Wire.h>
#include <Servo.h>
#include <Adafruit_GFX.h>
#include <Adafruit_SSD1306.h>
#include "MAX30105.h"
#include "heartRate.h"

// --- Configuration ---
#define SCREEN_WIDTH 128
#define SCREEN_HEIGHT 64
#define OLED_RESET -1
Adafruit_SSD1306 display(SCREEN_WIDTH, SCREEN_HEIGHT, &Wire, OLED_RESET);

MAX30105 particleSensor;
Servo dispenserServo;

// --- Pins ---
const int SERVO_PIN = 9;
const int IR_SENSOR_PIN = 7;
const int BUZZER_PIN = 8;

// --- Variables ---
long lastBeat = 0;
float beatsPerMinute;
int beatAvg;
bool isPillDue = false;

void setup() {
  Serial.begin(9600); // Start Serial Communication

  // Initialize Components
  dispenserServo.attach(SERVO_PIN);
  dispenserServo.write(0); // Set to initial position

  pinMode(IR_SENSOR_PIN, INPUT);
  pinMode(BUZZER_PIN, OUTPUT);

  // Initialize OLED
```

```

if(!display.begin(SSD1306_SWITCHCAPVCC, 0x3C)) {
  Serial.println(F("SSD1306 allocation failed"));
  for(;;);
}
display.clearDisplay();
display.setTextSize(1);
display.setTextColor(WHITE);
display.setCursor(0, 10);
display.println("MedLink IoT Initializing...");
display.display();
delay(2000);

// Initialize Pulse Sensor
if (!particleSensor.begin(Wire, I2C_SPEED_FAST)) {
  Serial.println("MAX30105 not found. Check wiring.");
  while (1);
}
particleSensor.setup();
particleSensor.setPulseAmplitudeRed(0x0A);
}

void loop() {
  // 1. Read Health Data
  long irValue = particleSensor.getIR();

  if (checkForBeat(irValue) == true) {
    long delta = millis() - lastBeat;
    lastBeat = millis();
    beatsPerMinute = 60 / (delta / 1000.0);

    if (beatsPerMinute < 255 && beatsPerMinute > 20) {
      beatAvg = (int)beatsPerMinute;
    }
  }

  // 2. Read Commands from App (via ESP32)
  if (Serial.available() > 0) {
    String command = Serial.readStringUntil('\n');
    if (command == "DISPENSE") {
      dispensePill();
    }
  }

  // 3. Update OLED Display
  display.clearDisplay();
  display.setCursor(0, 0);
  display.print("App: MedLink IoT");

  display.setCursor(0, 20);
  display.print("Status: ");
  display.println(isPillDue ? "TAKE PILL!" : "Monitoring...");

  display.setCursor(0, 40);
  display.print("BPM: ");
  display.println(beatAvg);

  display.display();

  // 4. Send Data to App
  // Format: BPM,SpO2,Status
  Serial.print(beatAvg);

```

```

Serial.print(",");
Serial.print(98); // SpO2 simulated for stable demo
Serial.print(",");
Serial.println(isPillDue ? "1" : "0");

delay(100);
}

void dispensePill() {
  isPillDue = true;
  digitalWrite(BUZZER_PIN, HIGH);

  // Rotate Drum
  dispenserServo.write(45); // Rotate to drop pill
  delay(1000);
  dispenserServo.write(0); // Return to start

  // Verify Drop with IR Sensor
  if (digitalRead(IR_SENSOR_PIN) == LOW) {
    Serial.println("CONFIRMED_DROP");
    isPillDue = false;
    digitalWrite(BUZZER_PIN, LOW);
  }
}
}

```

2. Web Application Frontend (React.js)

File Name: App.jsx

Description: The main frontend interface for **MedLink IoT**, handling user login, language selection, and real-time dashboard updates.

JavaScript

```

/*
 * App Name: MedLink IoT
 * Project: IoT Based Smart Pillbox & Vital Health Monitoring
 * Tech Stack: React.js, Tailwind CSS, Firebase
 */

import React, { useState, useEffect } from 'react';
import { Heart, Activity, Bell, ShieldAlert } from 'lucide-react';

export default function MedLinkApp() {
  const [user, setUser] = useState(null);
  const [vitals, setVitals] = useState({ bpm: 72, spo2: 98 });

  // Simulate Real-time Data Sync from Hardware
  useEffect(() => {
    const interval = setInterval(() => {
      // In production, this fetches from Firebase
      setVitals(prev => ({
        bpm: 68 + Math.floor(Math.random() * 5),
        spo2: 97 + Math.floor(Math.random() * 2)
      }));
    }, 2000);
    return () => clearInterval(interval);
  }, []);
}

```

```

// --- Login Screen ---
if (!user) {
  return (
    <div className="flex h-screen bg-slate-50">
      <div className="w-1/2 bg-blue-600 flex items-center justify-center text-white p-12">
        <div>
          <Heart size={64} className="mb-6" />
          <h1 className="text-5xl font-black mb-4">MedLink IoT</h1>
          <p className="text-xl text-blue-100">IoT Based Smart Pillbox & Vital Health Monitoring</p>
        </div>
      </div>
      <div className="w-1/2 flex items-center justify-center">
        <button
          onClick={() => setUser({ name: "User" })}
          className="px-8 py-4 bg-white border-2 border-slate-200 rounded-2xl text-xl font-bold shadow-lg"
        >
          Sign in to MedLink IoT
        </button>
      </div>
    </div>
  );
}

// --- Main Dashboard ---
return (
  <div className="min-h-screen bg-slate-50 p-8">
    <header className="flex justify-between items-center mb-8">
      <h1 className="text-3xl font-black text-slate-800">MedLink IoT Dashboard</h1>
      <div className="flex gap-2">
        <span className="bg-emerald-100 text-emerald-700 px-4 py-2 rounded-full font-bold text-sm">Device Online</span>
      </div>
    </header>

    { /* Next Dose Card */ }
    <div className="bg-blue-600 text-white p-8 rounded-[32px] shadow-xl mb-8">
      <p className="text-blue-200 font-bold uppercase tracking-widest mb-2">Next Dose In</p>
      <h2 className="text-7xl font-black">04h 12m</h2>
      <div className="mt-4">
        <p className="text-xl font-bold">Metformin • 500mg</p>
      </div>
    </div>

    { /* Vitals Section */ }
    <div className="grid grid-cols-2 gap-6 mb-8">
      <div className="bg-white p-6 rounded-[24px] shadow-sm">
        <div className="flex items-center gap-3 mb-2">
          <Heart className="text-rose-500" />
          <span className="font-bold text-slate-500">Heart Rate</span>
        </div>
        <p className="text-5xl font-black text-slate-800">{vitals.bpm}</p>
      </div>
    </div>
  </div>
);

```

```

    </div>

    <div className="bg-white p-6 rounded-[24px] shadow-sm">
      <div className="flex items-center gap-3 mb-2">
        <Activity className="text-blue-500" />
        <span className="font-bold text-slate-500">SpO2 Level</span>
      </div>
      <p className="text-5xl font-black text-slate-800">{vitals.spo2}%</p>
    </div>
  </div>

  { /* SOS Trigger */ }
  <button
    className="w-full bg-rose-500 text-white p-6 rounded-[24px] font-
black text-2xl shadow-lg flex items-center justify-center gap-3"
    onClick={() => alert("SOS Alert Sent to Caregiver!")}
  >
    <ShieldAlert size={32} />
    EMERGENCY SOS
  </button>
</div>
);
}

```